

# Keratinocyte Growth Supplement-defined (KGS-d) Catalog #2162

# **Product Description**

Keratinocyte Growth Supplement-defined (KGS-d) is a sterile, concentrated (100X) solution which contains growth factors, hormones, and proteins necessary for the culture of normal keratinocytes. The supplement is formulated (quantitatively and qualitatively) to provide an optimally balanced environment that selectively supports the growth of normal keratinocytes *in vitro*.

# **Product Use**

KGS-d is for research use only. It is not approved for human or animal use, or for application in *in vitro* diagnostic procedures.

### Storage

Store KGS-d at  $-20^{\circ}$ C; once added to medium, store at  $4^{\circ}$ C and protect from light. Do not refreeze after thawing.

# Shipping

Dry ice.

### **Instructions for Use**

For best results, use KGS-d in conjunction with Keratinocyte Medium-defined-b (KM-d-b, Cat #2111-b).

Thaw KGS-d at 37°C. Gently tilt the tube several times to ensure complete mixing. Spray the medium bottle and supplement tube with 70% ethanol and wipe to remove excess liquid. In a sterile field, remove the caps without touching the interior threads with fingers. Add KGS-d to the medium and mix well. Since several components are light-labile, the medium should not be exposed to light for extended periods. We do not recommend warming medium in a 37°C water bath prior to use. When stored in the dark at 4°C, the reconstituted medium is stable for one month.

*Caution: If handled improperly, some components of the medium may present a health hazard. Take appropriate precautions when handling it, including the wearing of protective clothing and eyewear. Dispose of properly.*