



**Human Hair Inner Root Sheath Cell Pellet  
(HHIRSC Cell Pellet)**  
Catalog #CP2430

Quantity: 5 million cells

**Product Description**

Human Hair Inner Root Sheath Cell Pellet (HHIRSC Cell Pellet) is prepared from early passage Human Hair Inner Root Sheath Cells. Each pellet contains 5 million cells and can be used for a variety of applications including PCR, western blotting, genomic DNA library construction, and gene expression profiling. Cell pellets from ScienCell Research Laboratories are convenient and cost effective for researchers as it eliminates the need to acquire expensive tissues.

**Applications**

- Genomic DNA
- cDNA
- Total RNA
- microRNA
- Protein Lysates

**Product Use**

HHIRSC Cell Pellet is for research use only. It is not approved for human or animal use, or for application in *in vitro* diagnostic procedures.

**Storage**

Store cell pellet immediately at -80°C upon arrival. Avoid freeze/thaw cycles.

**Shipping**

Dry ice.